

**QUICK DEFROSTING
EASY & CLEAN**

izidefrost+

DEFROSTING VACUUM

PROMOTE MAINTENANCE OF YOUR FREEZERS

-20°C / -80°C / -150°C

The essential tool for any laboratory with a fleet of freezers. Promote a more regular defrosting process, avoiding the need to turn off the freezer and transfer its contents. In doing so, extend the lifespan of your freezers and reduce your electricity consumption by **-30% à -50%**

● *Dimensions: L375mm x H53mm x W390mm - Weight: 9kg
Airflow: 4000l/m - Sound level: 64dba*

Professional device

Technical specifications specifically optimized for the aspiration of ice and snow at -80°C

Stainless steel scraper

Scrape and break the frost

25L tank

Vacuum the frost and ice

HEPA filter at the air outlet

Prevent any contamination in case of vacuuming a sample or sample residue





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.