

TECHNICAL FEATURE

7NG 35GBLP 00 905



The  :

- Comfortable and airy surface
- Without staples nor seams
- Abrasion-resistant
- Easily adjustable

ADJUSTABLE CHAIRS / POLYURETHANE

Chrome, 5-branch, scuff-resistant, reinforced base with articulated glides
Chrome gas-lift column with height adjustment from 410 to 540 mm and shock absorber
Ergonomic seat in black KANGO soft polyurethane
Articulated backrest in black soft polyurethane, adjustable in depth and height



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.