

## TECHNICAL FEATURE

7NN 10GHJR 02 905



### STOOLS / POLYURETHANE

Chrome, 5-branch, scuff-resistant, reinforced base with heavy-duty nylon casters  
Chrome gas-lift column with height adjustment from 580 to 820 mm and shock absorber  
Height adjustable polyamide foot-ring  
Round seat in black polyurethane

#### The :

- Without staples nor seams
- Abrasion-resistant
- Easily adjustable
- Polyvalent



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.