

# KERN 318-22

**KERN**

Set of milligram weights, E2, 1 mg - 500 mg, Platelet, stainless steel polished, in plastic case



Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	328-222-420
-------------------------------	-------------

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

## Measuring System

Test weight - nominal value	1 mg - 500 mg
OIML-compliant	✓
Weight OIML class	E2

## Construction

Dimension (W×D×H)	136×53×22 mm
Shape weight	Platelet
Material	stainless steel polished
Number of weights	13
Case of weight set	in plastic case

## Services (optional)

Article number for conformity assessment	952-350
Article number for DAkkS calibration	962-350
Article number for guaranteed plus-deviation	921-350

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	130×60×22 mm
Net weight	0,09 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.