

KERN ODC241



Digital microscopy brought up to date - tablet with integrated camera for optimal observation and digital documentation of the sample



Functions

Supported operating system: Android

Category

Brand: Optics

Packaging & shipping

Delivery: 24 h
Dimensions packaging (WxDxH): 275 x 230 x 85 mm
Gross weight: 1,050 kg
Net weight: 650 g
Shipping method: Parcel service

Construction

Colour/monochrome: Farbe
Dimensions (WxDxH): 238 x 51 x 206 mm
Resolution: 5 MP
Sensor: CMOS
Sensor size: 1/2.5"

Piktogramme

STANDARD





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.