

Test weights class F2

Class F2 Milligram weights, flat polygonal sheet, aluminium/German silver

Test weight material: Aluminium 1 mg – 5 mg/German silver 10 mg – 500 mg
 Container material: Lined plastic



Milligram weight			+		Container		+		DAkkS certificate		=	Package price
KERN		Tol ± mg			KERN				KERN			KERN
338-01	1 mg	0,06			347-009-400			962-451				
338-02	2 mg	0,06			347-009-400			962-452				
338-03	5 mg	0,06			347-009-400			962-453				
338-04	10 mg	0,08			347-009-400			962-454				
338-05	20 mg	0,10			347-009-400			962-455				
338-06	50 mg	0,12			347-009-400			962-456				
338-07	100 mg	0,16			347-009-400			962-457				
338-08	200 mg	0,20			347-009-400			962-458				
338-09	500 mg	0,25			347-009-400			962-459				

Class F2 Individual weights, cylindrical shape, finely turned stainless steel

Test weight material: finely turned stainless steel
 Container material: Lined plastic or wooden box (337-150-200 and 337-160-200)



Individual weight			+		Container		+		DAkkS certificate		=	Package price
KERN		Tol ± mg			KERN				KERN			KERN
337-01	1 g	0,3			347-030-400			962-431				
337-02	2 g	0,4			347-030-400			962-432				
337-03	5 g	0,5			347-030-400			962-433				
337-04	10 g	0,6			347-050-400			962-434				
337-05	20 g	0,8			347-050-400			962-435				
337-06	50 g	1,0			347-070-400			962-436				
337-07	100 g	1,6			347-070-400			962-437				
337-08	200 g	3,0			347-080-400			962-438				
337-09	500 g	8,0			347-090-400			962-439				
337-11	1 kg	16			347-110-400			962-441				
337-12	2 kg	30			347-120-400			962-442				
337-13	5 kg	80			347-130-400			962-443				
337-14	10 kg	160			347-140-400			962-444				
337-15	20 kg	300			337-150-200			962-445				
337-16	50 kg	800			337-160-200			962-446				

Alternative to plastic container:

Wooden boxes for individual weights. For more details on this, please see page 206.



Class F2 Test weights, stainless steel, stackable

Test weight material: finely turned stainless steel
 Box material: Wood



Test weight				+		Container		+		DAkkS certificate		=	Package price
KERN		Tol ± mg	Dim. Ø x H			KERN				KERN			KERN
337-141	10 kg	160	137x132 mm			337-141-200			962-444				
337-151	20 kg	300	137x217 mm			337-151-200			962-445				
337-161	50 kg	800	198x250 mm			337-161-200			962-446				



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

