

Test weights class F2

Class F2 Weight sets, cylindrical shape, finely turned stainless steel

Test weight material: 1 Individual weights - finely turned stainless steel,
milligram weights - aluminium 1 mg - 5 mg/German silver 10 mg - 500 mg
Case material: Lined plastic. Milligram weights 1 mg - 500 mg in plastic box.



Weight set		+	DAkkS certificate		=	Package price
KERN			KERN			KERN
338-22	1 mg - 500 mg		962-450			
333-024	1 mg - 50 g		962-401			
333-034	1 mg - 100 g		962-402			
333-044	1 mg - 200 g		962-403			
333-054	1 mg - 500 g		962-404			
333-064	1 mg - 1 kg		962-405			
333-074	1 mg - 2 kg		962-406			
333-084	1 mg - 5 kg		962-407			
334-024	1 g - 50 g		962-415			
334-034	1 g - 100 g		962-416			
334-044	1 g - 200 g		962-417			
334-054	1 g - 500 g		962-418			
334-064	1 g - 1 kg		962-419			
334-074	1 g - 2 kg		962-420			
334-084	1 g - 5 kg		962-421			

Class F2 Weight sets, cylindrical shape, finely turned stainless steel

Test weight material: Individual weights - finely turned stainless steel,
milligram weights - aluminium 1 mg - 5 mg/German silver 10 mg - 500 mg
Case material: Wood. Milligram weights 1 mg - 500 mg in plastic box.



Weight set		+	DAkkS certificate		=	Package price
KERN			KERN			KERN
338-22	1 mg - 500 mg		962-450			
333-02	1 mg - 50 g		962-401			
333-03	1 mg - 100 g		962-402			
333-04	1 mg - 200 g		962-403			
333-05	1 mg - 500 g		962-404			
333-06	1 mg - 1 kg		962-405			
333-07	1 mg - 2 kg		962-406			
333-08	1 mg - 5 kg		962-407			
333-09	1 mg - 10 kg		962-408			
334-02	1 g - 50 g		962-415			
334-03	1 g - 100 g		962-416			
334-04	1 g - 200 g		962-417			
334-05	1 g - 500 g		962-418			
334-06	1 g - 1 kg		962-419			
334-07	1 g - 2 kg		962-420			
334-08	1 g - 5 kg		962-421			
334-09	1 g - 10 kg		962-422			



You can create your own individual weight set yourself:

It contains only the weights which you need for testing purposes. KERN will customise your own personal box out of plastic, wood or aluminium. For more details on this, please see page 207.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.