

PRODUCT DATA SHEET

Standing: 2024-03-12

LAUDA Varioshake VS 20 OH

Shaker 230 V; 50/60 Hz

Part Number: L003131

Features

- Shaker with overhead rotation, heavy duty, ideal for mixing powders and pasty samples
- Speed display
- Electronic speed control, stepless, gentle start-up
- Clear control panel for easy operation with rotary knobs
- Capacity at uniform loading: 12 vessels, round or square
- Vessels up to 110 mm Ø and 270 mm height, height adjustable clamp holders
- Drive by an overload-protected DC motor
- Rotating rack with ball bearings on both sides with safe drive via slide hub
- Height-adjustable holding devices made of stainless steel for safe, firm clamping of the vessels
- Constant speed, independent of load
- Excellent for continuous use
- Outer housing made of electrolytically galvanized sheet steel, powder-coated



Technical Features

Ambient temperature range	10 ... 40 °C
Power consumption max.	0.1 kW
Shaking frequency	1 ... 20 1/min
Load max.	20 kg
Overall dimensions (WxDxH)	770 x 700 x 715 mm
Weight	68 kg
Power plug	Power cord with angled plug (CEE7/7)
Power supply	230 V; 50/60 Hz

Reserve technical changes



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.