




COUNTERTOP LABORATORY FRIDGE SOLID DOOR 47L



- ✓ **+2° C to +10° C Operating Temperature Range**
Conforms to MHRA guidelines
- ✓ **Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside unit
- ✓ **Audible Warning Alarms**
Alerts user should internal temperature deviate from set parameters

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444411348	LSFSR47DC-UK	Countertop Laboratory Fridge Solid Door 47L White



FULL SPECIFICATION – LSFSR47DC-UK COUNTERTOP LABORATORY FRIDGE SOLID DOOR 47L WHITE

Description	Countertop
SKU Code	444411348
EAN Code	5052263113480
GDPA Code	LSFSR47DC-UK
Gross Capacity (L)	47
Net Capacity (L)	30
Unit External Width (mm)	475
Unit External Depth (mm)	445
Unit External Height (mm)	565
Product Weight (kg)	17.5
Packed Width (mm)	490
Packed Depth (mm)	530
Packed Height (mm)	620
Packed Weight (kg)	21.5
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	409 x 354 x 425
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	+2° C to +10° C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	905
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet	Yes
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	2
Number of Shelf Positions	3
Shelf Size	416 x 200
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Automatic
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	No
Insulation Thickness (mm)	40
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

