



Solid Door



Glass Door

## Countertop Pharmacy Essential Refrigerator 107L

GDPA code: PESR107UK (Solid), PEGR107UK (Glass)

SKU code: 444410629 (Solid), 444410630 (Glass)

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 2 year warranty as standard



The Experts in  
Medical Refrigeration

## Specification: Countertop Pharmacy Essential Refrigerator 107L

<b>Gross Capacity (litres)</b>	107
<b>Net Capacity (litres)</b>	103
<b>Unit External Width (mm)</b>	495
<b>Unit External Depth (mm)</b>	620
<b>Unit External Height (mm)</b>	830
<b>Product Weight (kg)</b>	31.5kg (solid), 36.5kg (glass)
<b>Packed Width (mm)</b>	560
<b>Packed Depth (mm)</b>	850
<b>Packed Height (mm)</b>	655
<b>Packed Weight (kg)</b>	35kg (solid), 40kg (glass)
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	403mm W x 700mm H x 400mm D
<b>Temperature Range</b>	2°C to 8°C
<b>Colour</b>	White/green
<b>External Finish</b>	Painted steel
<b>Internal Finish</b>	HIPS*
<b>Type of plastic</b>	SRL800 (HIPS/ABS mix)
<b>Door Finish</b>	Painted steel
<b>Door Type</b>	Solid (solid), Glass (glass)
<b>Door Style</b>	Hinged reversible
<b>Door Open Depth (mm)</b>	1050mm
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	Yes
<b>Lock</b>	Yes, 2 Keys
<b>Adjustable Feet (number &amp; position)</b>	1 adjustable front foot
<b>Porthole</b>	No
<b>Porthole Size</b>	N/A
<b>Porthole Location</b>	N/A
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	Yes
<b>Remote Alarm Contacts</b>	Yes
<b>Audible Alarms</b>	Yes
<b>Visual Alarms</b>	Yes
<b>Number of Shelves</b>	2
<b>Number of Shelf Positions</b>	4
<b>Shelf Size (mm) (WxDxH)</b>	427mm W x 6mm H x 374mm D
<b>Shelf Material</b>	Wire
<b>Number of Baskets</b>	1
<b>Basket Material</b>	Wire
<b>Basket Size (mm) (WxDxH)</b>	—
<b>Interior Light</b>	No
<b>Light Type</b>	N/A
<b>Light Activation</b>	N/A
<b>Fan</b>	Yes
<b>Spark Free</b>	No
<b>Auto or Manual Defrost</b>	Auto
<b>Controller Type</b>	Digital, external
<b>Number of Probes</b>	1
<b>Probe Type</b>	Air
<b>Climate Class</b>	N
<b>Refrigerant Type</b>	R600a
<b>CFC Free</b>	Yes
<b>HFC Free</b>	No

**Specification:** Countertop Pharmacy Essential Refrigerator 107L

<b>HCFC Free</b>	Yes
<b>Ammonia Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Noise Level (dB)</b>	<43dB
<b>Energy Consumption (kwh)</b>	—
<b>Voltage Supply</b>	230
<b>Aux Output</b>	Yes
<b>Insulation Thickness (mm)</b>	35mm
<b>Insulation Type</b>	Cyclopentane
<b>Blowing Agent</b>	Yes
<b>Anti Microbial Handle</b>	No
<b>Optional Drawer System Available</b>	No
<b>Wall Brackets</b>	No
<b>Digital Lock</b>	No



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.