



## Countertop Pharmacy Essential Refrigerator 47L

GDPA code: PESR47DC-UK

SKU code: 444411286

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms




The Experts in  
Medical Refrigeration



## Specification: Countertop Pharmacy Essential Refrigerator 47L

Gross Capacity (litres)	47
Net Capacity (litres)	30
Unit External Width (mm)	480
Unit External Depth (mm)	505
Unit External Height (mm)	565
Product Weight (kg)	17.5
Packed Width (mm)	490
Packed Depth (mm)	530
Packed Height (mm)	620
Packed Weight (kg)	21
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	400 x 230 x 415
Temperature Range	+2°C to +8°C (set point tolerance + or - 2)
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	905
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	4 - All adjustable
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	No
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	1
Number of Shelf Positions	3
Shelf Size (mm) (WxDxH)	416 X 200
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Basket Size (mm) (WxDxH)	396 x 165
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	Yes - Axial
Spark Free	No
Auto or Manual Defrost	None
Controller Type	Digital, External
Number of Probes	1
Probe Type	Air
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes

**Specification:** Countertop Pharmacy Essential Refrigerator 47L



<b>HCFC Free</b>	Yes
<b>Ammonia Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	230v
<b>Insulation Thickness (mm)</b>	40
<b>Insulation Type</b>	High Density Polyurethane Foam
<b>Blowing Agent</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes
<b>Optional Drawer System Available</b>	No
<b>Wall Brackets</b>	Yes
<b>Digital Lock</b>	No





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.