



## Large Pharmacy Refrigerator Solid Door 1200L PSR1200UK

PSR1200UK

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display showing air and core temperatures
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ Temperature monitoring software records air and load temperatures constantly and also holds a full 7 day alarm history
- ✓ USB data download
- ✓ Power failure alarm

## SPECIFICATION

Dimensions	1300mm W x 1995mm H x 840mm D
Weight	200kg
Product name	Large Pharmacy Refrigerator Solid Door 1200L PSR1200UK
GDPA code	PSR1200UK
SKU code	UK 444440368, EU 444442594
External dimensions	1300mm W x 1995mm H x 840mm D
Internal dimensions	1220mm W x 1386mm H x 653mm D
Packed dimensions	1480mm W x 2195mm H x 990mm D
Product weight	200kg
Packed weight	220kg
Gross capacity	1200L
External finish	Stainless steel

Colour	White/green
Internal finish	Stainless steel
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door finish	Painted stainless steel
Door style	Hinged
Open door depth	1145mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	Fitted castors
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	10
Shelf material	Wire
Number of shelf positions	40
Shelf size	645mm W x 523mm D
Interior light	N
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	ST
Refrigerant type	R134a
CFC free	Y
HFC free	Y
HCFC free	Y
Ammonia free	Y
UK 3 pin plug	Y
Voltage supply	230

Spark free	N
Auxiliary output	Y
Insulation thickness	60mm
Insulation type	Cyclopentane
Battery low alarm	Y (visual)
Air/load/defrost sensor fail alarm	Y (visual)
Compressor alarm	Y (visual)
Door open alarm	Y (visual and audio)
Air tempertaure alarm	Y (visual and audio)
Mains failure alarm	Y (visual and audio)
Antimicrobial handle	N
Optional drawer system available	N



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.