



LABCOLD™

medical & scientific refrigeration



19 litre Portable Vaccine Carrier

This 19 litre Labcold portable vaccine and sample carrier is designed to solve the problem of cold chain compliance on the move as it can be powered by both mains and vehicular power sources.

There is an easy to read digital temperature display and because it is accurate to one decimal place, the vaccine carrier can be calibrated so you can be confident that the temperature displayed is as experienced by the contents of the carrier.

Using a thermoelectric, solid state cooling system there are no dangerous liquids that could escape or moving parts that could get damaged.

❄ Supplied with two power leads. 230 volt for the mains and the other designed to fit 12v power sockets in vehicles

❄ Precise digital controller with real time and set point temperature display means the carrier can be validated in accordance with DH Green Book chapter 3 requirements

❄ Smooth, easy clean plastic exterior, powder coated interior

❄ Compact and lightweight



The carrier comes complete with a convenient shoulder strap and the interior features a divider to keep the contents in place during transit.

Purpose designed for any professional who needs to give vaccines in the field, this carrier can be used in conjunction with Labcold data loggers for extra temperature security.

Please be aware this product is not designed to be run for extended periods of time nor to hold temperature when not connected to power source.

Part Number	RPDF0019F
Temperature Range	+2°C to +10°C
Exterior Dimensions (HxWxD)	355 x 465 x 320mm
Weight	6.5kg
Energy Consumption	1.5 kW/24hr
Carry Strap	✓
Digital temperature control	✓
Temperature display	✓
Moulded feet for stability in transit	✓
12 volt DC and 240 AC volt operation	✓



12v and mains power points



Clear digital temperature display



Hygienic easy clean exterior





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.