

RPFR21043 Solid Door Advanced Pharmacy Refrigerator

CLAIM
YOUR
FREE
2
YEAR
WARRANTY
LABCOLD



Robustly constructed from steel with a stainless steel interior this refrigerator from the **Labcold** Advanced Pharmacy Range is designed to provide close temperature control for temperature sensitive medicines in the pharmacy.

In addition, remote alarm contacts can be fitted as well as access ports. Contact us for bespoke racking inventory systems or specialised applications.



The air is continually circulated internally to ensure that where ever product is placed inside the cabinet, it stays within range. There is no internal 'step' and the interior is hygienic stainless steel making it easy to clean and stock. The 4 heavy duty shelves are fully adjustable and more shelves can be purchased as an optional extra. The robust doors are independently lockable and for increased security and convenience are self closing.

The cabinet is fitted with heavy duty, lockable castors making it easy to manoeuvre and clean behind which means this Advanced Pharmacy Refrigerator is ideal for demanding environments.

| Part Number | RPFR21043 |
|-----------------------------|------------------------------------|
| Capacity (Litres) | 650 |
| Style | Freestanding |
| Exterior Dimensions (HxWxD) | 1985 x 710 x 800mm |
| Energy Consumption | 2.43kWh/24h |
| Operating Temperature | +5°C |
| Temperature Set Point Range | +2°C to +8°C |
| Interior Construction | Stainless Steel |
| Exterior Construction | White polyester coated sheet steel |
| Shelves | 4 |
| Maximum shelf loading kg | 70 |
| Shelf dimensions (W x D) | 500 x 620 mm |
| Door Hinging | Right side hinging only |
| Door Style | Solid |
| Internal Light | ✓ |
| Lockable Castors | ✓ |
| Door Lock | ✓ |
| Access Port | Optional Extra |
| Digital Temperature Control | ✓ |
| Digital Temperature Display | ✓ |
| High Temperature Alarm | ✓ |
| Low Temperature Alarm | ✓ |
| Door Open Alarm | ✓ |
| Maximum Temperature Display | ✓ |
| Minimum Temperature Display | ✓ |
| Auto Defrost | ✓ |
| Fan Assisted Cooling | ✓ |
| Low GWP Refrigerant | ✓ |
| 2 Year Warranty (UK Only) | ✓ |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.