



# RAFR05204 Advanced Under Counter Freezer

Part Number	RAFR05204
Capacity (Litres)	115
Style	Under Counter
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Energy Consumption	1.5kWh/24h
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -22°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W x D)	460 x 4075mm 460 x 240mm
Door Hinging	Right (can be reversed at time of order)
Self Closing Door	✓
Door Lock	✓
Levelling Feet	4
Access Port	✓
Remote Alarm Contacts	✓
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Manual defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
2 Year Warranty*	✓

Robustly constructed with a stainless Steel exterior this freezer from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products .



Insulated with Polyurethane Water Base Insulation the refrigerant is also ammonia, HCFC and CFC free and uses a GWP friendly refrigerant.

Fitted with a sturdy door lock, this Advanced freezer is ideal for demanding environments such as busy laboratories and features a manual defrost putting you in control.

Remote alarm contacts and an access port are also fitted as standard.

\*Parts only warranty outside of the UK





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.