



# LABCOLD™

medical & scientific refrigeration



## RLPL04044 Basic Laboratory Refrigerator

Constructed from enamelled steel, this refrigerator from the Labcold Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the Labcold brand and support.

The fridge has two non adjustable shelves and a polycarbonate drawer at the bottom of the refrigerator. The adjustable feet ensures the fridge can be stabilised and with auto defrost and an internal light, the refrigerator has some useful convenience features.

To put you in control, the temperature can be easily adjusted with a dial located on the right hand side of the interior. Labcold will by default set the set point to +4°C.

Please note, this refrigerator is not sparkfree and can not have a lock or digital temperature display fitted to it.

This Basic Refrigerator is ideal for economic cold storage of items that are not volatile and where precise storage temperature is not product critical. As you would expect from Labcold, this Basic Refrigerator comes complete with a full one years parts and labour warranty.\*



Part Number	RLPL04044
Capacity (Litres)	91
Style	Underbench
Exterior Dimensions (HxWxD)	842 x 475 x 448mm
Weight	22kg
Energy Consumption	0.29 kWh/24h
Noise Level	39 dB(A)
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +10°C
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Shelves	2 + 1 bin
Sparkfree	✗
Door Hinging	Right (can be reversed at time of ordering)
Door Style	Solid
Feet	2 adjustable feet at front
Internal Light	✓
Auto Defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
1 year UK Parts and Labour Warranty	✓

\* Labour UK only





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.