



LABCOLD™

medical & scientific refrigeration



RLCG0600 Refrigerated Cooled Incubator

Part Number	RLCG0600
Capacity (Litres)	570
Style	Freestanding
Exterior Dimensions (HxWxD)	1980 x 720 x 840mm
Interior Dimensions (HxWxD)	1405 x 600 x 570mm
Weight	140kg
Internal Light	LED
Factory Set Point	+20°C
Operating Temperature	+5°C to +40°C
Interior Construction	Stainless Steel
Exterior Construction	Enamelled Steel
Door Style	Glass
Door Hinging	Right
Inner Doors	2
Shelves	4
Shelf dimension (WxDxH)	580 x 530 x 15mm
Spare shelves available	✓
Temperature display	✓
Microprocessor controlled temperature regulation	✓
Fan assisted temperature control	✓
Access Port	✓
Auto Defrost	X
Eco Friendly Refrigerant	✓

This 570 litre cooled incubator is ideal for microbiology or any application where precise and consistent temperatures at or below room temperature are required.

Digital
Microprocessor
Controller



With a robust stainless steel interior, enamelled steel exterior and fully adjustable shelves, this incubator is versatile and built to last.

Labcold Cooled Incubators are also available in 88 litre, 150 litre and 600 litre capacity

Most models have the choice of glass or solid door





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.