




COUNTERTOP LABORATORY FREEZER SOLID DOOR 39L



-18° C to -25° C Operating Temperature Range



Conforms to MHRA guidelines

Unique Anti Microbial Handle



Helps prevent the growth of bacteria and minimise the potential for cross contamination

Fitted Lock as Standard



Maintains security and ensures the unit can only be accessed by approved personnel

Wire Shelf



Enables optimum air flow inside unit

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444411346

LSFSF39DC-UK

Countertop Laboratory Freezer Solid Door 39L



LSFSF39DC-UK COUNTERTOP LABORATORY FREEZER SOLID DOOR 39L

Description	Countertop
SKU Code	444411346
EAN Code	5052263113466
GDPA Code	LSFSF39DC-UK
Gross Capacity (L)	39
Net Capacity (L)	35
Unit External Width (mm)	475
Unit External Depth (mm)	445
Unit External Height (mm)	565
Product Weight (kg)	18
Packed Width (mm)	490
Packed Depth (mm)	445
Packed Height (mm)	620
Packed Weight (kg)	21
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	370 X 345 X 404
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	-18° C to -25° C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	905
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet	Yes
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	1
Number of Shelf Positions	Fixed
Shelf Size	365 x 190
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Manual
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	No
Insulation Thickness (mm)	30
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

