




## High Accuracy External Thermistor Wifi Temperature Sensor Data Logger

GDPA code: EL-WIFI-DTP+GD Wifi Logger

SKU code: 444410733

- ✓ Thermistor probe temperature measurement range  
-40 to +125°C (-40 to +257°F)
- ✓ Wirelessly stream and view data via WiFi on PC or  
Lec Medical Cloud
- ✓ Easy sensor set-up using free PC software application
- ✓ View and analyse multiple sensors, including  
graphing of historic data
- ✓ Configurable high and low alarms with indicator
- ✓ Sensor memory stores all data even if WiFi is  
temporarily disconnected
- ✓ Cloud access - <https://lecmedical.wifisensorcloud.com>

**Specification:** High Accuracy External Thermistor Wifi Temperature Sensor Data Logger



<b>Dimensions</b>	82mm W x 70mm H x 36mm D
<b>Weight</b>	N/A
<b>Measurement range</b>	-40°C to 125°C
<b>Temperature display on panel</b>	Y
<b>Accuracy</b>	+/-0.6°C
<b>Number of readings shown</b>	2
<b>Number of probes</b>	2



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

