

# Pharmacy Refrigeration

## Portable Cooler Range

### Portable Vaccine Cooler 18L

GDPA code **VCP18UK 18L Portable Cooler**  
 SKU code **444410883**



### Key features

- ✓ +2°C to +8°C operating temperature range
- ✓ External LED temperature display
- ✓ Holds temperature inside between 2°C & 8°C for up to 30 hours
- ✓ Aluminium heat preserving foil lining
- ✓ Typical vaccine storage 12-13L
- ✓ Covered bung non slip feet
- ✓ Webbing adjustable shoulder strap
- ✓ Comes with 12 gel packs

### Specification: Portable Vaccine Cooler 18L

Height (mm)	310
Width (mm)	440
Depth (mm)	285
Weight (kg)	4.03
Gross Capacity (L)	18
Net Capacity (L)	8
Temperature Range	2 to 8 degrees celsius
Ambient Temperature Range	-10 to +40 degrees celsius
Temperature Hold Range (hours)	Upto 30
Aluminium Lined Sleeve	Yes
Handle Type	Plastic
Additional Shoulder Strap	Yes
Gel Packs Included	Yes
Number of Gel Packs	12
CE Certified	Yes
External Temperature Display	Yes
Internal Air Temperature Probe	Yes
Battery Included	Yes
Lockable	No



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

