

## Heating Stirrer Mantle (250ml)

**Product Code:** ES-302N

**Brand:** Medline Scientific

This mantle features an integral magnetic stirrer, which is independently controlled. The stirrer speed may be varied from 100 to 1000rpm and the temperature is regulated by an infinity bimetallic thermostat. Retort rod holder integrated into base of unit to accommodate our range of accessories.

### Specifications:

<b>Volume (ml)</b>	250
<b>Case Construction</b>	Steel, with chemical resistant coating
<b>Thermal Insulation</b>	Moulded silica powder
<b>Heating Element</b>	Nichrome wire
<b>Mantle</b>	Woven glass fibre mat
<b>Element Temperature</b>	Up to 450 °C
<b>Speed Range</b>	10 - 1000rpm
<b>Power Supply</b>	220V AC
<b>Watts</b>	150





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

