

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Miele ProCare Lab 10 MA

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

Washing and cleaning products (including solvent based products)

1.3. Details of the supplier of the safety data sheet

Address:

Miele & Cie. KG
Carl-Miele-Straße 29
D-33332 Gütersloh
Telephone no. +49 5241 89 0
Fax no. +49 5241 89 2090
www.miele-professional.com

E-mail address of person responsible for this SDS:

sida@drweigert.de

1.4. Emergency telephone number

GBK/ Infotrac: (USA domestic) 1 800 535 5053 or international +1 352 323 3500

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification (Regulation (EC) No. 1272/2008)

Classification (Regulation (EC) No. 1272/2008)

| | |
|---------------|------|
| Skin Corr. 1B | H314 |
| Eye Dam. 1 | H318 |
| Met. Corr. 1 | H290 |

2.2. Label elements

Labelling according to regulation (EC) No 1272/2008

Hazard pictograms



Signal word

Danger

Hazard statements

| | |
|------|--|
| H290 | May be corrosive to metals. |
| H314 | Causes severe skin burns and eye damage. |

Precautionary statements

| | |
|----------------|---|
| P280 | Wear protective gloves/protective clothing/eye protection/face protection. |
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. |

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER or doctor.
Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

Hazardous component(s) to be indicated on label (Regulation (EC) No. 1272/2008)

contains disodium metasilicate;dipotassium metasilicate

2.3. Other hazards

No special hazards have to be mentioned.

SECTION 3: Composition/information on ingredients

Hazardous ingredients

disodium metasilicate

CAS No. 6834-92-0
EINECS no. 229-912-9
Registration no. 01-2119449811-37-0004
Concentration $\geq 10 < 25$ %
Classification Xi, R37
C, R34

Classification (Regulation (EC) No. 1272/2008)
Skin Corr. 1B H314
STOT SE 3 H335

dipotassium metasilicate

CAS No. 10006-28-7
EINECS no. 233-001-1
Concentration $\geq 10 < 25$ %
Classification Xi, R37
C, R34

Classification (Regulation (EC) No. 1272/2008)
Skin Corr. 1B H314
STOT SE 3 H335

potassium hydroxide

CAS No. 1310-58-3
EINECS no. 215-181-3
Registration no. 01-2119487136-33
Concentration < 1 %
Classification Xn, R22
C, R35

Classification (Regulation (EC) No. 1272/2008)
Acute Tox. 4 H302
Skin Corr. 1A H314

Concentration limits (Regulation (EC) No. 1272/2008)
Skin Corr. 1A H314 ≥ 5
Skin Irrit. 2 H315 $\leq 0,5 < 2$
Skin Corr. 1B H314 $\leq 2 < 5$
Eye Irrit. 2 H319 $\leq 0,5 < 2$

SECTION 4: First aid measures

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

4.1. Description of first aid measures

General information

Remove contaminated, soaked clothing immediately and dispose of safely. Clean body thoroughly (bath, shower). In any case show the physician the Safety Data Sheet.

After inhalation

Ensure supply of fresh air. When spray fog inhaled, seek medical aid.

After skin contact

After contact with skin, wash immediately with plenty of water. Take medical treatment.

After eye contact

In case of contact with the eyes, rinse immediately for at least 15 minutes with plenty of water. Summon a doctor immediately.

After ingestion

If swallowed, seek medical advice immediately and show this container or label. Rinse mouth thoroughly with water. Let plenty of water be drunk in small gulps. Do not induce vomiting.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed

Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Extinguishing measures to suit surroundings

Non Suitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

5.3. Advice for firefighters

Special protective equipment for fire-fighting

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus.

Other information

Collect contaminated fire-fighting water separately, must not be discharged into the drains. Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing. Refer to protective measures listed in sections 7 and 8.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/groundwater.

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

6.3. Methods and material for containment and cleaning up

Pick up with absorbent material. Dispose of absorbed material in accordance with the regulations.

6.4. Reference to other sections

Refer to protective measures listed in sections 7 and 8.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid formation of aerosols. Observe the usual precautions for handling chemicals. Keep container tightly closed.

Advice on protection against fire and explosion

The product is not combustible.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep in original packaging, tightly closed. Storage rooms must be properly ventilated. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Storage class according to TRGS 510

Storage class according to TRGS 510 8B Non-combustible corrosive hazardous substances

7.3. Specific end use(s)

no data

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

There are not known any further control parameters.

8.2. Exposure controls

General protective and hygiene measures

Hold eye wash fountain available. Hold emergency shower available. Do not inhale gases/vapours/aerosols. Avoid contact with skin and eyes. Do not eat, drink or smoke during work time. Wash hands before breaks and after work. Clean skin thoroughly after work; apply skin cream.

Respiratory protection

If workplace limits are exceeded, a respiratory protection approved for this particular job must be worn.

Hand protection

Chemical resistant gloves

| | |
|----------------------|--------------|
| Appropriate Material | neoprene |
| Appropriate Material | butyl rubber |
| Appropriate Material | nitrile |

Eye protection

Safety glasses with side protection shield

Body protection

Clothing as usual in the chemical industry. Protective shoes

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

| | |
|---|------------------------|
| Form | liquid |
| Colour | colourless |
| Odour | characteristic |
| Odour threshold | |
| Remarks | not determined |
| pH value | |
| Value | appr. 14 |
| Melting point | |
| Remarks | not determined |
| Freezing point | |
| Remarks | not determined |
| Initial boiling point and boiling range | |
| Remarks | not determined |
| Flash point | |
| Remarks | Not applicable |
| Evaporation rate (ether = 1) : | |
| Remarks | not determined |
| Flammability (solid, gas) | |
| evaluation | not determined |
| Upper/lower flammability or explosive limits | |
| Remarks | not determined |
| Vapour pressure | |
| Remarks | not determined |
| Vapour density | |
| Remarks | not determined |
| Density | |
| Value | 1,48 g/cm ³ |
| Solubility in water | |
| Remarks | not determined |
| Solubility(ies) | |
| Remarks | not determined |
| Partition coefficient: n-octanol/water | |
| Remarks | not determined |
| Ignition temperature | |
| Remarks | not determined |
| Decomposition temperature | |
| Remarks | not determined |
| Viscosity | |
| Remarks | not determined |
| Explosive properties | |
| evaluation | not determined |
| Oxidising properties | |
| evaluation | None known |

9.2. Other information

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 / GB

Date revised: 26.11.2015

Print date: 26.11.15

Other information

None known

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

No hazardous reactions known.

Decomposition temperature

Remarks not determined

10.5. Incompatible materials

Strong exothermic reaction with acids. Corrodes aluminium.

10.6. Hazardous decomposition products

Irritant gases/vapours

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute oral toxicity

| | | | |
|---------|--|--|-------|
| Species | rat | | |
| LD50 | > 2000 | | mg/kg |
| Method | calculated value (Regulation (EC) No. 1272/2008) | | |

Acute oral toxicity (Components)

disodium metasilicate

| | | | |
|---------|------|--|-------|
| Species | rat | | |
| LD50 | 1150 | | mg/kg |

Acute dermal toxicity

Remarks not determined

Acute dermal toxicity (Components)

disodium metasilicate

| | | | |
|---------|--------|--|-------|
| Species | rat | | |
| LD50 | > 5000 | | mg/kg |

Acute inhalational toxicity

Remarks not determined

Skin corrosion/irritation

evaluation corrosive

Serious eye damage/irritation

evaluation corrosive

Sensitization

Remarks not determined

Subacute, subchronic, chronic toxicity

Remarks not determined

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 / GB

Date revised: 26.11.2015

Print date: 26.11.15

Mutagenicity

Remarks not determined

Reproductive toxicity

Remarks not determined

Carcinogenicity

Remarks not determined

Specific Target Organ Toxicity (STOT)

Remarks not determined

Aspiration hazard

No special hazards have to be mentioned.

Experience in practice

Inhalation may lead to irritation of the respiratory tract.

Other information

There is no data available on the product apart from the information given in this subsection.

SECTION 12: Ecological information

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

disodium metasilicate

| | | | |
|----------------------|---------------|---|------|
| Species | mosquito fish | | |
| LC50 | 2320 | | mg/l |
| Duration of exposure | 96 | h | |

disodium metasilicate

| | | | |
|----------------------|--------------------------------|---|------|
| Species | zebra fish (Brachydanio rerio) | | |
| LC50 | 210 | | mg/l |
| Duration of exposure | 96 | h | |

Bacteria toxicity (Components)

disodium metasilicate

| | | | |
|----------------------|------------------|---|------|
| Species | activated sludge | | |
| EC50 | > 100 | | mg/l |
| Duration of exposure | 3 | h | |

12.2. Persistence and degradability

General information

not determined

12.3. Bioaccumulative potential

General information

not determined

Partition coefficient: n-octanol/water

Remarks not determined

12.4. Mobility in soil

General information

not determined

12.5. Results of PBT and vPvB assessment

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

General information

not determined

12.6. Other adverse effects

General information

not determined

General information / ecology

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations for the product

Allocation of a waste code number, according to the European Waste Catalogue (EWC), should be carried out in agreement with the regional waste disposal company.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

Land transport ADR/RID

14.1. UN number

UN 1719

14.2. UN proper shipping name

CAUSTIC ALKALI LIQUID, N.O.S. (disodium metasilicate, dipotassium metasilicate)

14.3. Transport hazard class(es)

Class 8

Label 8

14.4. Packing group

Packing group III

Limited Quantity 5 I

Transport category 3

Tunnel restriction code E

Marine transport IMDG/GGVSee

14.1. UN number

UN 1719

14.2. UN proper shipping name

CAUSTIC ALKALI LIQUID, N.O.S. (disodium metasilicate, dipotassium metasilicate)

14.3. Transport hazard class(es)

Class 8

14.4. Packing group

Packing group III

14.5. Environmental hazards

no

IMDG-Code segregation group 18 Alkalis

Air transport ICAO/IATA

14.1. UN number

UN 1719

14.2. UN proper shipping name

CAUSTIC ALKALI LIQUID, N.O.S. (disodium metasilicate, dipotassium metasilicate)

14.3. Transport hazard class(es)

Miele ProCare Lab 10 MA

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 26.11.2015

Print date: 26.11.15

Class 8

14.4. Packing group

Packing group III

Information for all modes of transport

14.6. Special precautions for user

See Sections 6 to 8

Other information

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Ingredients (Regulation (EC) No 648/2004)

15 % or over but less than 30 %:

phosphates

Water Hazard Class (Germany)

Water Hazard Class WGK 1
(Germany)

Remarks Classification according to Annex 4 VwVwS

VOC

VOC (EU) 0 %

Other information

The product does not contain substances of very high concern (SVHC).

15.2. Chemical safety assessment

For this mixture a chemical safety assessment has not been carried out.

SECTION 16: Other information

R-phrases listed in Chapter 3

22 Harmful if swallowed.
34 Causes burns.
35 Causes severe burns.
37 Irritating to respiratory system.

Hazard statements listed in Chapter 3

H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.
H335 May cause respiratory irritation.

CLP categories listed in Chapter 3

Acute Tox. 4 Acute toxicity, Category 4
Skin Corr. 1A Skin corrosion, Category 1A
Skin Corr. 1B Skin corrosion, Category 1B
STOT SE 3 Specific target organ toxicity - single exposure, Category 3

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.