



## MB-800 Series Moisture Determination Balances

- Three-year industry leading warranty
- Readability as precise as 0.1 mg
- Touch screen display
- Automatic weighing chamber door



# Cole-Parmer® MB-800 Series Moisture Determination Balances

The MB-800 series of moisture balances combines state of the art features, an intuitive touch screen display, premium build quality, and Cole-Parmer's industry leading three-year warranty into an economical package that is an essential addition to any laboratory.

Onboard database storage and processing, and unparalleled connectivity in the form of RS-232, USB-A, USB-B, and Wi-Fi make the MB-800 series moisture determination balances the perfect choice for all your moisture analysis needs.

## Technical Specifications

Specification	MB-800-50.H	MB-800-50	MB-800-210
Max Capacity	50 g	50 g	210 g
Readability	0.1 mg	1 mg	1 mg
Calibration Type	External	External	External
Weigh Pan Dimensions	ø90 mm	ø90 mm	ø90 mm
Heater Type	Infrared	Infrared	Infrared
Max Temperature	160 °C	160 °C	160 °C
Device Dimensions	13 x 8 x 7.8 in (338 x 206 x 197 mm)		
Connectivity	2x RS-232, 2x USB-A (interchangeable), USB-B, Wi-Fi		

## Ordering Information

Description	Series No.	Model No.	Ordering No.
50 g, 0.1 mg Moisture Determination Balance with Touch Screen, High Resolution	MB-800	MB-800-50.H	10100-42
50 g, 1 mg Moisture Determination Balance with Touch Screen	MB-800	MB-800-50	10100-43
210 g, 1 mg Moisture Determination Balance with Touch Screen	MB-800	MB-800-210	10100-44



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.