

# Circulaire® Powder Containment



PCB 1800

**BRITISH  
MANUFACTURED**



**FIVE YEAR  
WARRANTY**



**GUARANTEED  
PROTECTION**



**MAXIMUM  
OPERATOR  
SAFETY**



**The Circulaire® Powder Containment Booth offers first class operator protection from both Powders and Particulate.**

The unit is designed specifically for control of medium to heavy powders whilst dispensing from containment drums to a balance or smaller container. The high velocity rear extraction offers excellent containment by utilising H14 HEPA Filters, 99.997% efficient at 0.3µm and filtering air up to 1250m<sup>3</sup>/hr, achieving an air velocity of up to 0.7m/s.

## Circulaire® Powder Weighing

**The Circulaire® Powder Weighing Cabinet utilises Main & Outlet ULPA Filters to meet COSHH compliance. The unit is specifically designed for control of medium to heavy powders whilst dispensing to a balance or smaller container.**

The Granite work surface ensures accurate weighing for sensitive applications and the removable and hinged baffles make for easy cleaning.



PWC  
1100

### Size Specification (W x D x H)

PCB 1800	PWC 1100
1800mm x 750mm x 2400mm	1100mm x 700mm x 1284mm



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.