


CAT. NO.....:	450-0020/-0045/-0080	
PORE SIZE (µm).....:	0.2/0.45/0.8	
UPPER CAP. (mL).....:	500	
RECEIVER CAP. (mL)::	500	
MEMBRANE COLOR.....:	WHITE	
MEMBRANE GRID COLOR (-0020/-0045/-0080)::	NONE/GREEN/BLACK	
MEMBRANE DIA (mm)::	75	
	INCHES	MM
A (UPPER I.D.).....:	3.80	96.5
B HEIGHT.....:	9.8	249
C DIA (BASE).....:	3.86	98.0
D OVERALL.....:	4.0	102
NOM. WT. (gms).....:	155	

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	GREEN PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HOPE
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	CN
VENT PLUG	CA

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS 00-0450-XXXX-1	REV.
CHK'D	PVC		6/25/99		0699

DESCRIPTION:

STERILIZATION FILTER UNIT



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.