


| | | | | | | |
|----------------------|----------|-------|----------|------|----------|------|
| CAT. NO..... | 565-0020 | | 568-0020 | | 566-0020 | |
| UPPER CAP. (mL).... | 150 | | 250 | | 500 | |
| RECEIVER CAP. (mL).. | 150 | | 250 | | 500 | |
| PORE SIZE (µm)..... | 0.2 | | 0.2 | | 0.2 | |
| MEMBRANE COLOR..... | WHITE | | WHITE | | WHITE | |
| MEMBRANE DIA (mm).. | 50 | | 50 | | 75 | |
| | INCHES | MM | INCHES | MM | INCHES | MM |
| A (UPPER I.D.)..... | 3.15 | 80.0 | 3.15 | 80.0 | 3.80 | 96.5 |
| B HEIGHT..... | 6.15 | 156.2 | 7.9 | 200 | 9.8 | 249 |
| C DIA (BASE)..... | 3.50 | 88.9 | 3.50 | 88.9 | 3.86 | 98.0 |
| D OVERALL..... | 3.7 | 94 | 3.7 | 94 | 4.0 | 102 |
| NOM. WT. (gms)..... | 78 | | 95 | | 155 | |

| MATERIALS | |
|-----------------|-------------|
| UPPER CONTAINER | PS |
| COVER | PS |
| COLLAR | BLUE PS |
| RECEIVER | PS |
| STORAGE CLOSURE | WHITE HDPE |
| ADAPTOR | WHITE LDPE |
| GASKET | FOAMED LDPE |
| MEMBRANE | αPES |
| VENT PLUG | CA |

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

| | | | | | | |
|-----------------------------------|-----|-------|----------|--------------------|------|------|
| MADE | PB | APP'D | JMC | PSS 00-056X-0020-1 | REV. | 1107 |
| CK'D | JMC | | 11/29/07 | | | |
| DESCRIPTION: αPES FILTER UNITS | | | | | | |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.