



### Unique Features:

- Secure and simple setup, with a factory paired Bluetooth adapter leading to plug-and-play integration.
- With its statistics and totalization modes and clock, the SF40A/BT ensures accurate and compliant reporting. These features complement OHAUS weighing instruments' application modes and provide increased functionality for laboratory, industrial, and educational applications.
- With the ability to print in 13 different languages, SF40A/BT can print in the language of the connected device to support universal operation.



# SF40A/BT Bluetooth Impact Printer

Model	SF40A/BT
Item number	30944359
Number of Columns	40 columns
Print Speed	Approximately 1,2 lines per second
Application Modes	Totalization, Statistics (Statistics printout includes number of samples, average value, standard deviation, minimum, maximum, min-max differential, and total sample weight sum)
Paper Roll (W x Dia.)	Standard paper 57,5 mm x Ø 60 mm max.
Power Requirements	AC adapter Input: 100 – 240 VAC, 50 – 60 Hz AC adapter output: 12 VDC, 2,5 A
Interface	1 x RS232, 2 x Bluetooth Adaptors
Ambient Temperature	0 °C – 40 °C
Printer Dimensions (W x L x H)	140 x 205 x 75 mm
Shipping Dimensions (W x L x H)	255 mm x 205 mm x 150 mm
Wireless standard*	2,0
Active Communication Range*	10 m (without obstacles)
Net Weight	824 g (paper roll included)
Shipping Weight	0,8 kg

## Other Standard Features and Equipment

Status indicator, 40 characters per line, perform zero and tare operations directly from the printer. Full functionality with Explorer™, Navigator™, MB Series Moisture Analyzers and 'T' Series indicators. Statistics and Totalization not supported on Scout™ Series. Multi-Language Printing: English, Spanish, German, Italian, French, Chinese, Korean, Japanese, Russian, Polish, Portuguese, Hungarian and Czech.


## Compliance

- **SF40A printer:** Product Safety: IEC/EN 61010-1, : UL61010-1; CAN/CSA C22.2 61010-1
- **Electromagnetic Compatibility:** FCC Part 15 Class A; Canada ICES-003 Class A; IEC/EN 61326-1 Class B, basic environments
- **Compliance Marks:** CSA, CE
- **Bluetooth adapter certification:** Anatel (for Brazil), CE RED (for EU), FCC (for US), IC (for Canada), iDA (for Singapore), KCC (for Korea), MoC (for Israel), NBTC (for Thailand), NCC (for Taiwan China), NOM (for Mexico), QAS SIRIM (for Malaysia), RCM (for Australia and New Zealand), SRRC (for China), Telec (for Japan), Bluetooth SIG.

## Accessories

- Paper Roll x 2 (included) ..... 12120799
- Paper Rotary Axis (included) ..... 30063920
- Ink Ribbon (included) ..... 30529322
- Cable Adapter (included)..... 30059316
- Cable Mini USB ..... 30944363

## Outline Dimensions





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.