




BAC 130

Information

OUTER DIMENSIONS	800 x 600 x H 715 mm
INNER DIMENSIONS	620 x 410 x H 520 mm
INNER VOLUME	132 Litres
K COEFFICIENT	0,36W/m2.K
CURB WEIGHT	31 kg
INSULATION	Polyurethane foam
WALLS	Food-graded polyethylene
OUTER SIZE	Half EUR pallet : 600x800mm - Euronorm plastic crate: 600x400mm




Specifications



Chest and lid made of food-graded monobloc rotomolded polyethylene



8 molded handles



Inner convection groovings



Embedded label holder




Eutectic plate molded shelf



Molded bossage enabling stacking of BAC



Chest/lid molded hinge, 115° opening



Monobloc frame silicone seal



Embedded closing latch - sealable with a padlock/seal

Accessories



Eutectic plate TOP 130



Dry ice drawer



Removable castor base fitted with 4 wheels Ø 100 mm: 2 fixed and 2 swivel




Removable castor base fitted with 4 wheels Ø 100 mm: 2 fixed and 2 swivel with brakes



Removable castor base fitted with 4 swivel wheels Ø 100mm : 2 swivel and 2 swivel with brakes



Castor base can be fixed on request



Several colors available under conditions on request



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.