



* Optical fiber (CL-11.1) NOT included

Model:
CLD-01

Typology:
ILLUMINATOR

Description:
LED cold light generator

Light source	High power 14 W LED.
Chromaticity coordinates (x,y)	0.320, 0.320 (pure white, 6,300 K).
Optical power (with double-arm)	Illuminance (at 10 cm distance): 70,000 lux.
Power supply	Brightness control potentiometer. External power supply provided. Input: 110/240 Vac 50/60 Hz; output: 6 Vdc 5 A.
Available optical fibers	- CL-11.1: Double-arm optical fiber guide. Each arm fitted with focusable lens adapter. Arm length 500 mm. - CL-17.1: Polarizing filters (pair) for CL-11.1. - CL-12: Ring optical fiber guide: Guide length 700 mm. Guide diameter 16 mm. Diameter of the fixing ring: 55 mm. The circular end is suitable for all stereo series.
Dimensions	Length: 200 mm Width: 120 mm Weight: 4.1 Kg



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

