



Model:

CL-41

Typology:

ILLUMINATOR

Description:

Double flexible arm (L = 60 cm) with double X-LED 3W.

Chromaticity coordinates (x,y):

0.320, 0.320 (pure white)

Optical power:

Luminous flux: 400 lumen

Illuminance (at 10 cm distance): 170000 lux

Aperture angle:

25°

Power supply:

Brightness control potentiometer.

External power supply provided.

Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

