



96-LED ring polarized light illuminator,
multi-plug

Chromaticity and LED	Coordinates (x,y)	(0.32, 0.32)
	Color Temperature (K)	6,300
	LED number	96
	Brightness control	Manual
	Lifetime (hours)	> 65,000
	Max. required power (W)	6
Optical power	Luminous flux (lumen)	480
	Illuminance (at 10 cm distance, lux)	9,000 (on 6 cm diameter)
Power Supply for Illumination	Type	External
	Connector	Jack, 1.4 mm
	Power plug type	Multi-plug (EU, UK, US)
	Input voltage	100/240 Vac, 50/60 Hz
	Output voltage	12 Vdc 1.5 A
Accessories Included	Fixing ring	Yes
	Rotating Polarizing filter	Yes
	Analyzer filter	Yes



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

