



*Binocular stereozoom microscope head,
6.7x...45x*

| | | |
|---|------------------------------------|--|
| Observation Method - Transmitted Light | Brightfield | As optional |
| | Darkfield | As optional |
| | Simple polarized light | As optional |
| Observation Method - Incident Light | Brightfield | Yes |
| | Simple polarized light | As optional |
| Head | Type | Binocular |
| | Construction material | Plastic mold / Aluminum gears |
| | Inclination | 45° |
| | 360° rotating | Yes |
| | Interpupillary distance (mm) | 51-75 |
| | Fixing screw for eyepieces | Yes |
| | Head holder internal diameter (mm) | 76 |
| Tube inner diameter (mm) | 30 | |
| Eyepieces | Field number (mm) | 23 |
| | Magnification | 10x |
| | Planar type | Yes |
| | Micrometric scale | As optional |
| | Diameter of micrometer glass (mm) | 26 |
| | High eyepoint (for glass wearers) | Yes |
| | Dioptric adjustment | Yes |
| | Rubber cup | Yes |
| Retractable protections | Yes | |
| Objective | Optical system | Greenough |
| | Anti-fungus treatment | Yes |
| | Zoom Type | Parfocal achromatic |
| | Working distance (mm) | 110 |
| | Standard magnifications | 6.7x-45x |
| | Zoom ratio | 6.72 |
| | Zoom click stops | At 0.67x, 0.8x, 1x, 1.2x, 1.5x, 2x, 2.5x, 3x, 3.5x, 4x, 4.5x |
| Accessories Included | Dust cover | Yes |
| | User Manual | Digital version (downloadable) |
| Product Dimensions | Head height (mm) | 240 |
| | Head width (mm) | 167 |
| | Head depth (mm) | 200 |
| Product Weight | (kg) | 1.8 |



Coarse focusing system

| | | |
|-----------------------------|------------------------------------|--------------------------------|
| Focusing System | Type | Coarse |
| | Coarse total travel (mm) | 50 |
| | Adjustable tension | Yes |
| | Head holder internal diameter (mm) | 76 |
| Accessories Included | User Manual | Digital version (downloadable) |
| | | |
| Product Dimensions | Total height (mm) | 100 |
| | Total width (mm) | 150 |
| | Total depth (mm) | 250 |
| Product Weight | (kg) | 1 |



Simple stand, pillar type

| | | |
|--|--|--------------------------------|
| Observation Method - Incident Light | Brightfield | As optional |
| | Simple polarized light | As optional |
| Stand | Type | Pillar |
| | Pillar diameter (mm) | 32 |
| Stage | Black/white contrast plate for incident ill. | Yes, diameter 95 mm |
| | Specimen clips | Yes |
| Accessories Included | Dust cover | Yes |
| | User Manual | Digital version (downloadable) |
| Product Dimensions | Total height (mm) | 275 |
| | Base height (mm) | 25 |
| | Base width (mm) | 210 |
| | Base depth (mm) | 270 |
| Product Weight | (kg) | 1.4 |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.