



Model:  
**CL-14**

Typology:  
**ILLUMINATOR**

Description:  
*LED ring light illuminator.*

<b>Cromaticity coordinates (x,y):</b>	0.320, 0.320 (pure white), 6,300 K.
<b>Optical power:</b>	Luminous flux: 240 lumen. Illuminance (at 10 cm distance): 8,000 lux (on 3 cm diameter).
<b>Power supply:</b>	Brightness adjustment knob. Main power switch 0/1. External power supply provided. Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1,500 mA.
<b>Suitable models:</b>	LAB, SZM and SZN.
<b>Accessories:</b>	Fixing ring for stereo heads.



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

