



144-LED ring illuminator, EU plug

<b>Chromaticity and LED</b>	Coordinates (x,y)	(0.32, 0.32)
	Color Temperature (K)	6,300
	LED number	144
	Brightness control	Manual
	Quadrant ON/OFF control	Yes
	Lifetime (hours)	> 65,000
	Max. required power (W)	2
<b>Optical power</b>	Luminous flux (lumen)	120
	Illuminance (at 10 cm distance, lux)	6,000 (on 10 cm diameter)
<b>Power Supply for Illumination</b>	Type	Internal
	Power plug type	EU plug
	Input voltage	100/240 Vac, 50/60 Hz, 150mA
	Output voltage	10-12Vdc
<b>Accessories Included</b>	Fixing ring	Yes



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

