

Tuttnauer

LEADING IN
STERILIZATION
SINCE
1925

T-CLASSIC

9-Litre Portable Autoclave

Making your process easy, quick and affordable.

The T-Classic One, is a lightweight, portable steam autoclave designed as a cost-effective sterilization solution. Fully automatic with an easy one touch operation, delivering short cycles.

Ideally designed to meet portability needs of a beauty clinic or academic facility.

The T-Classic One can also be used in combination with another autoclave or as a back up.

Tuttnauer Performance with the simplicity and affordability your small practice needs



T-Classic One

Instruments sterilized and ready quickly



The T-Classic One offers competitively quick unwrapped cycle times, in as short as 18 minutes*. The fast point-of-use sterilization makes your turnaround a priority, increasing efficiency, while eliminating bottlenecks.

*Cycle times may vary, depending on load & weight.



The T-Classic One is a simple one step operated device. With one touch, the T-Classic One will begin an automatic cycle until completion. An electronic detector will turn off the power if water level is insufficient. The perfect autoclave for limited space.



The T-Classic One is the lightest tabletop autoclave available, weighing only 5.35 kg allowing for easy portability, while requiring minimal counter space.

T-CLASSIC



The **T-Classic One** offers multiple benefits for users. The Autoclave has 9 litre capacity, with a 126°C sterilizing temperature. The internal furniture consists of a basket as standard, however can also be supplied with cassette holders for instruments.

Technical Specifications

Maximum Width Overall	340mm
Maximum Height Overall	335mm
Maximum net weight	5.35 kg
Chamber Diameter/Height	210mm/230mm
Chamber capacity	9 Litres
Maximum Instrument length	228mm
Maximum Load weight	3kg
Sterilisation Temperature	126 °C
Operating Pressure	1.6 Bar
Minimum Sterilisation time	11 minutes
Normal Cycle Time	22 minutes
Voltage	230 Volts
Wattage	1500 Watts
Frequency	50 – 60 Hertz
Furniture	Basket



Top view, easy access



Basket for easy loading and unloading



Easy to use control panel

Model	Inner Chamber Dimensions ø x D (mm)	Volume (Litres)	Shipping Weight (kg)	Shipping Dimensions WxDxH (mm)	Single Phase Power Supply
T-Classic One	210 x 230	9	8.5	400 x 400 x 620	220-240V, 50-60Hz 1500W, 10A

LTR330-0128-A



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.