

T-CLASSIC

Podiaclave



The T-Classic Podiaclave has a 9 litre capacity, with a 134 °C sterilizing temperature. The internal furniture consists of a basket as standard and the unit also comes complete with a printer. The portability of T-Classic Podiaclave means that it is ideal for those practitioners who domiciliary work and need to take a sterilizer perform with them.



Model Number	T-Classic Podiaclave
Maximum Width Overall	340mm
Maximum Height Overall	335mm
Maximum net weight	6.40 kg
Chamber Diameter/Height	210mm/230mm
Chamber capacity	9 Litres
Maximum Instrument length	228mm
Maximum Load weight	3kg
Sterilisation Temperature	134 °C
Operating Pressure	2.2 Bar
Minimum Sterilisation time	3 minutes
Cycle Time Normal	20 minutes
Voltage	230 Volts
Wattage	1500 Watts
Frequency	50 – 60 Hertz
Printer	Yes
Temperature and Pressure Gauges	N/A
Furniture	Basket



Top view, easy access



Basket for easy loading and unloading



Easy to use control panel

Model	Inner Chamber Dimensions ø x D (mm)	Volume (Litres)	Shipping Weight (kg)	Shipping Dimensions WxDxH (mm)	Single Phase Power Supply
T-Classic Podiaclave	210 x 330	9	9.25	400 x 400 x 620	220-240V, 50-60Hz 1500W, 10A



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.