

Classic Portable Autoclave



 **Prestige Medical**®
Integrated decontamination solutions



Prestige Medical Classic – with more than 1 million sold in 80 countries worldwide, you can rely on our reputation for quality and performance.

Technical Specification

Classic Portable Autoclave							
Model Number	210001	210002	210004	210047	210048	210134	
Max width - overall	340mm	340mm	340mm	340mm	340mm	340mm	340mm
Max height - overall	335mm	420mm	420mm	335mm	420mm	420mm	335mm
Max net weight	4.0kg	4.8kg	4.8kg	4.0kg	4.8kg	4.8kg	5.4kg
Chamber diameter	210mm	210mm	210mm	210mm	210mm	210mm	210mm
Chamber height	236mm	328mm	328mm	236mm	328mm	328mm	236mm
Chamber capacity	9 litres	12 litres	12 litres	9 litres	12 litres	12 litres	9 litres
Max instrument length	228mm	290mm	290mm	228mm	290mm	290mm	228mm
Max load weight	3kg	4kg	4kg	3kg	4kg	4kg	3kg
Sterilizing temperature	126°C	126°C	126°C	121°C	121°C	121°C	134°C
Min sterilizing time	11 minutes	11 minutes	11 minutes	28 minutes	28 minutes	28 minutes	3 minutes
Cycle time - normal	22 minutes	22 minutes	22 minutes	70 minutes	70 minutes	70 minutes	20 minutes
Operating pressure	1.4 bar	1.4 bar	1.4 bar	1.05 bar	1.05 bar	1.05 bar	2.2 bar
Voltage	230v	230v	230v	230v	230v	230v	230v
Wattage	1250w	1250w	1250w	1250w	1250w	1250w	1250w
Frequency	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz
Printer							
Temperature & Pressure Gauges			✓				✓
Furniture	219708	219720	219720	219708	219720	219720	219708

All Classic units are supplied complete with electrical cord set and instruction manual as appropriate.



Cassette x 3 &
Carrier - ref: 219706

Cassette
ref: 219291

Cassette
Carrier
ref: 219740



Cord Sets
UK ref: 219258
Euro ref: 219297
UL ref: 219299

A wide choice of accessories - allows you to specify the model to suit your exact needs.

Standard
Basket
ref: 219708



Reliable and trusted worldwide

Designed to give you the most cost effective sterilization solution with a wide choice of models and accessories that are guaranteed to meet your needs.

Steam sterilization, or autoclaving, is the preferred method of eliminating the risk of cross infection through contaminated instruments. It is the only method recommended throughout the world by healthcare professionals including the World Health Organisation (WHO) and the Centre for Disease Control (CDC).

Sterilizing Indicator Strips



126°C ref: 219277
121°C/134°C ref: 259277

Extended
Basket
ref: 219720



Lifter Fork
ref: 219294



Green Gaskets x 1
ref: 219500





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.