

## FLOOR CUPBOARDS



CZ709046ZKXX



CZ703530ZKXX



CZ189046ZKXX

### Floor Cupboards - Dimensions/Codes

Size H.W.D.	Shelves	Galvanised Code	Price	Extra Shelves	Code	Price
1800.1200.460	3	CZ181246ZKXX	£508.51	W.1200 x D.460	CUPDSHLF1246	£52.16
1800.900.460	3	CZ189046ZKXX	£431.47	W.900 x D.460	CUPDSHLF9046	£45.04
1200.900.460	2	CZ129046ZKXX	£366.27	W.900 x D.460	CUPDSHLF9046	£45.04
900.900.460	1	CZ909046ZKXX	£290.41	W.900 x D.460	CUPDSHLF9046	£45.04
900.460.460	1	CZ904646ZKXX	£226.40	W.460 x D.460	CUPDSHLF4646	£34.38
700.900.460	1	CZ709046ZKXX	£265.52	W.900 x D.460	CUPDSHLF9046	£45.04
700.350.300	1	CZ703530ZKXX	£192.03	W.350 x D.300	CUPDSHLF3530	£28.45

## MOBILE CUPBOARDS



### Mobile Cupboards - Dimensions/Codes

Size H.W.D.	Shelves	Code	Price
1040.900.460	1	CZ109046MKXX	£461.10
840.900.460	1	CZ849046MKXX	£423.17

## FLOOR CHESTS



### Floor Chests - Dimensions/Codes

Size H.W.D.	Code	Price
610.1170.460	CZ611146ZKXX	£490.73
510.610.340	CZ516134ZKXX	£312.93

## WALL CUPBOARDS



### Wall Cupboards - Dimensions/Codes

Size H.W.D.	Shelves	Code	Price
570.850.255	1	CZ855725ZKXX	£251.29

## STANDS



### Stands - Dimensions/Codes

Size H.W.D.	Code	Price
543.900.460	CZST9046ZKXX	£116.16
543.460.460	CZST4646ZKXX	£103.13
543.350.300	CZST3530ZKXX	£99.57



It is vital that employers and employees comply with the control of hazardous substances in all working environments. This range of cupboards and cabinets enables the separate storage and transportation of hazardous substances as demanded by the CoSHH regulations.

- “No snag” handles with 2-point locking.
- Adjustable spill-retaining galvanised shelves for superior spill management\*.
- New regulation labels enclosed to be fitted by customers to suit regulations.
- Integral liquid-tight sump encourages the correct handling of spills through the use of absorbents\*.
- Wall cupboards are pre-drilled for mounting (fixings not included).
- Floor chests have a welded integrated sump, angled lid to discourage article placement and a hasp and staple lock - padlock not supplied.
- Stands raise cupboards by 543mm with an adjustable foot for uneven surfaces.
- Order quantities may apply, please speak to the sales team.
- Powder coated with Germ Guard Active Technology Anti -Bacterial Paint.
- Reaction to Fire Classification EN 13501-1.
- Cupboard shelves have a UDL of 50kg.

\*Except wall mounted cupboards.





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.