


# Q55 Sonicator

**The Q55 is a compact and cost effective ultrasonic processor that will occupy less bench space than any unit on the market.**

This model is effective for standard cell disruption and many other small volume applications. Probes are available in three different sizes.

## FEATURES:

- Smallest unit available
- Thumb-switch or continuous operation
- Simple and effective operation



Stand sold separately.

## PART NO. Q55 INCLUDES:

- Generator
- Converter
- 1/8" (3mm) dia. probe
- Power cable
- Converter cable
- Wrench set

## TECHNICAL SPECIFICATIONS:

<b>Power Rating:</b>	55 watts
<b>Frequency:</b>	20kHz
<b>Dimensions (W x L x H):</b>	8 x 7.5 x 5.75 in. (203 x 190 x 146 mm)
<b>Voltage:</b>	110V, 50/60 Hz

Specify desired voltage for export.

## Optional Accessories

Part #	Description
460	Support Stand with Converter Holder
432A	Sound Enclosure with Converter Holder Exterior Dimensions (W x H x D) 20 x 12 x 12 in. (508 x 254 x 254 mm)
CL188	Replacement Converter



Probes	Part # 4423	Part # 4422	Part # 4435
Processing Volume*	200µl - 5ml	500µl - 15ml	10ml - 50ml
Tip Diameter	5/64" (2mm)	1/8" (3mm)	1/4" (6mm)
Amplitude (µm)	200	180	120



#CL188



#460



#432A



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.