

INTEGRA Biosciences PIPETBOY Pro

INTEGRA

Highly accurate flow control is achieved simply by applying variable finger pressure on the pipettor trigger buttons, changing between drop-by-drop dispensing or fast motorised pipetting at speeds of up to 10mL/s.

Supplied with a spare filter, two batteries, mains adapter and wall mount.

- Integrated LED illumination
- Fine control of pipetting via finger pressure enables accurate drop-by-drop dispensing
- Quiet and strong pump for high pipetting speed of up to 10mL/s
- Lightweight and ergonomic design for maximum comfort
- Removable and autoclavable pipette holder module

Code	Colour
I156400	White
I156401	Pink
I156402	Green
I156403	Orange
I156404	Blue

Accessories

Code	Description
I156450	Charging station for parking and charging the PIPETBOY pro
I155065	Stand for PIPETBOY and VACUBOY
I156632	Adapter UK, 240V AC
I151020	Pipette mount for holding the pipette in the nose piece, silicone
I156151	Battery set, 2 AAA batteries, rechargeable



SLS Lab Basics Pipette Filler

SLS LAB BASICS

Lightweight, made from UV resistant material, features adjustable valves. Filling and dispensing speed is easily set using the same hand by pressing the large mode button. The bright LED's will indicate which of the three modes has been selected, High, Low or Gravity Mode.

5 coloured nose cones are included. User labelling area provides customisation, instant recognition and/or task designation. Additional features include hand neutral design, concave padded triggers, a 1mL pipette adapter to increase stability of small pipettes, a long life Li-Ion battery that allows 20 hours of continuous usage, low battery indicator and two holders as standard.

- For 1-100mL glass or plastic pipettes
- Uses standard 0.45µm hydrophobic filters
- Complete with universal voltage and plug faced charger
- Comes with bench stand and magnetic/fixable wall mount

Code	Description
SLS7010	Pipette filler
PIP6042	Sterile filter 0.45µm



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.