

# TeleScoop Pro

The Extendable Sampler For Liquid Sampling

Volume 250ml to 2000ml

**New  
sizes**

The all new TeleScoop Pro is perfect for sampling liquids and suspensions.

The lightweight aluminium telescopic rods can be easily extended to the ideal length. A wide range of fittings can be quickly connected to ensure effortless sampling.

All Attachments are interchangeable with all Telescopic Rods. Our innovative Quick-Click Change System allows attachments to be changed in just a few seconds - no tools are required.

### Ordering:

1. Choose your Telescopic Rod
2. Choose your attachments



## 1. Telescopic Rods

- Wide range of sizes
- New improved robust design
- Quickly extends to the required length in seconds
- Lightweight



Part No.	Description	Length	Number of Sections	Material
5065A-1001	Telescopic Rod 1000mm	1000mm	1	Aluminium & plastic
5065A-1201	Telescopic Rod 1200mm	750 to 1250mm	2	Aluminium & plastic
5065A-2501	Telescopic Rod 2500mm	1370 to 2500mm	2	Aluminium & plastic
5065A-2801	Telescopic Rod 2800mm	1050 to 2800mm	4	Aluminium & plastic
5065A-3001	Telescopic Rod 3000mm	1235 to 3000mm	3	Aluminium & plastic
5065A-4501	Telescopic Rod 4500mm	1720 to 4500mm	3	Aluminium & plastic
5065A-6001	Telescopic Rod 6000mm	1845 to 6000mm	4	Aluminium & plastic

Liquid Samplers



**Quick-Click  
Change  
System**

No Tools Required

## 2. TeleScoop Pro Attachments



### Angular Beaker 5065B

The beaker can be easily adjusted to the required angle



### Wrist Strap 3111A

Ensures the sampler is not lost if it is accidentally dropped during use



### Pendulum Beaker 5065C

Always swings to the upright position



### Stainless Pendulum Beaker 5065C

Sampling beaker made from 316L stainless steel. Always swings to the upright position



### Bottle Holder 5065D

Suitable for bottles up to 95mm diameter. Price includes one 1000ml HDPE bottle & lid



### Stainless Bottle Holder 5065D

Sampling bottle made from 316L stainless steel

Bottle supplied with lid and clamp

Part No.	Description	Nominal Volume	Material
5065B-252	Angular Beaker	250ml	PP
5065B-502	Angular Beaker	500ml	PP
5065B-1002	Angular Beaker	1000ml	PP
5065B-1001B	Angular Beaker	1000ml	316L Stainless Steel
5065B-2002	Angular Beaker	2000ml	PP
5065C-252	Pendulum Beaker	250ml	PP
5065C-502	Pendulum Beaker	500ml	PP
5065C-1002	Pendulum Beaker	1000ml	PP
5065C-1001B	Pendulum Beaker	1000ml	316L Stainless Steel
5065C-2002	Pendulum Beaker	2000ml	PP
5065D-502	Bottle Holder	500ml	HDPE/PP Bottle
5065D-503	Bottle Holder	500ml	316L Stainless Steel
5065D-1002	Bottle Holder	1000ml	HDPE/PP Bottle
SPARES			
5065G-1B	Elbow Fitting	n/a	PBT
8126H-500	Spare Bottle	500ml	HDPE body with PP lid. Box of 48
8126H-1000	Spare Bottle	1000ml	HDPE body with PP lid. Box of 24
A836-500	Spare Bottle	500ml	316L Stainless Steel
3111A-01	Wrist Strap	n/a	PET

### Elbow Fitting 5065G-1B

The multi-angle Elbow Fitting allows the operator to achieve the perfect angle for sampling and is included with all of the attachments shown above (except the wrist strap).





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.