

	<h1>Product Specification</h1>	Specification Number	8121P-8
		Revision	01

Part No. 8121P-8

Description Narrow Mouth Bottle - 8ml

Material of Construction: PP (polypropylene)

- Translucent
- FDA and USP Class VI compliant




Method of Construction: Cap - Single piece, injection moulded  
Body - Single piece, injection moulded

Overall Height: 45mm

Maximum Diameter: 25mm

Nominal Weight: 6g

BSE/TSE Free: Yes

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 25-Jan-13
Page 1 of 1		



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

