


**1. Part Number:** A417-5

**2. Photograph:**



**3. Product Information:**

Description	Stainless Container 5L with External GL45 Thread
Materials of Construction:	<ul style="list-style-type: none"> <li>- Body: 316L stainless steel</li> <li>- Lid: 316L stainless steel</li> <li>- Clamp: 304 stainless steel</li> <li>- Gasket: Silicone (FDA acceptable grade)</li> <li>- Cap: 316L stainless steel</li> <li>- Cap Seal: Silicone (FDA acceptable grade)</li> </ul>
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	5 Litres
Brimful Volume:	5.7 Litres
Dimensions	A = 119mm B = 98mm C = 316mm D = 164mm
Weight:	2.7 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Note:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Oct-20



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

