

1. Part Number: A426-25

2. Photograph:



3. Product Information:

Description:	Mini-Pot 25ml
Method of Construction:	Crevice free
Material of Construction:	316L Stainless Steel
Nominal Volume:	25ml
Brimful Volume:	29ml
Overall Height (body & lid):	36mm
Body Height:	33mm
Maximum Body Diameter:	38mm
Lid Diameter:	41mm
Nominal Weight:	110g
Finish:	Bright Polish (Ra < 0.8 microns)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: If you require container with a watertight seal please see our containers type A836

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Mar-21



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

