



	<h1>Product Specification</h1>	Specification Number	A566-10
		Revision	02

Part No. A566-10
Description TC Container 10L
Nominal Volume: 10L
Material of Construction: Body: 316L stainless steel
Side Handles: 304 stainless steel
Lid: 316L stainless steel
Clamp: 316L stainless steel (may be 304 stainless steel)
Gasket: Silicone (conforms to FDA CFR 177.2600)



Method of Construction: Welds ground & polished, crevice free interior
Surface Finish: Better than 0.5 microns Ra
Dimensions: A = 119mm
B = 98mm
C = 366mm
D = 240mm
E = 210mm
Handles: Two side handles welded on to patch plates
Graduation Marks: None
Pressure Rating: None
Nominal Weight (body, lid, clamp & gasket): 3.65 Kg

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 29-Aug-19
----------------------------	---	---



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

