

	<h1>Product Specification</h1>	Specification Number	8209H-20S
		Revision	01

Part No. 8209H-20S

Description SteriWare Long Handled Spoon (Straight)

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004



Method of Construction: Single piece, injection moulded

Bowl Volume: 20ml

Handle Length: 285mm

Bowl Length: 88mm

Bowl Width: 38mm

Nominal Weight: 22g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)


Method of Sterilisation: Gamma Irradiation (min dose 25kGy)

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 13-Jul-18
Page 1 of 1		



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.