

1. Part Number: 8404H-1000S

2. Photograph:



3. Product Information:

Description	Wide Mouth HDPE Bottle 1000ml - Sterile
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Packing Environment:	Class 10,000 Cleanroom
Nominal Volume:	1000ml
Brimful Volume:	1075ml
Nominal Height:	201mm
Diameter:	92mm
Internal Neck Diameter:	52mm
Weight:	100g
Individually Bagged:	Yes (in a PE/PET laminate film bag)
Quantity per box:	30
Method of Sterilisation:	Gamma Irradiation (25-45 kGy)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Jul-21





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

