

1. Part Number: A201-150

2. Photograph:



3. Product Information:

Description:	Pallet Knife with Wooden Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Rosewood
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	146mm
Blade Width:	23mm
Overall Length:	252mm
Nominal Weight:	39g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jan-22



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

