

**1. Part Number:** A296-1

**2. Photograph:**



**3. Product Information:**

Description	Stainless Steel Bailer 1L
Materials of Construction:	316L stainless steel
Nominal Volume:	1 Litre
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Handle Length:	115mm
Max. Bowl Diameter:	170mm
Bowl Height:	70mm
Nominal Weight:	435g
Graduation Marks:	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Dec-19



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.