

1. Part Numbers: A686-15

2. Photograph:



3. Product Information:

Description	Bucket with Spout 15L
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free interior Lugs for handle are spot welded to exterior of bucket
Surface Finish:	Better than 1.0 microns Ra
Nominal Volume:	15 L
Brimful Volume:	15.2 L
Max. width over spout:	355mm
Internal Diameter at Top:	310mm
Diameter at Base:	202mm
Bucket Height:	333mm
Graduations:	No
Skirt Band:	Yes
Nominal Weight:	3.5kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

